

i-Pro™ Personal Respiratory Protective Mask – Fitting guide

Preparations for use



1. Wash hands before use and ensure the mask is undamaged and free of any obstructions or foreign bodies.



2. Cup the mask in your hand, avoiding contact with the filter, allowing the straps to hang freely below your hand.

Fitting the mask



3. Position the mask under your chin with the nosepiece up. Pull the bottom strap over your head resting it around the neck below the ears.



4. Pull the top strap over your head and position it above the ears, at the back of your head.



5. Hold the front of the mask to the face and pull the top elastic loop at the side to tighten the mask – take care to not contact the filter media while putting on the mask.



6. Again, holding the mask to the face, pull the bottom loop elastic at the side to tighten the mask.



7. To check the mask to face seal, hold two fingers above the mask to check no air is escaping around the seal as you exhale.



8. Then hold both hands underneath the mask to check no air is escaping around the seal as you exhale. Be careful not to disturb the position of the mask. If air leaks, adjust position of straps and make certain mask edges fit snugly against the face.