

# Intersurgical Guedel Airways

## New ISO 5364:2016 Specified Colours

### Colour change requirement

In order to comply with the colour coding system specified in ISO 5364:2016, it was necessary for Intersurgical to update four sizes of our Guedels. The sizes that have changed are as follows:

- ISO size 3.5; from Pink to Light Green
- ISO size 6.5; from White to Brown
- ISO size 7.0; from Yellow to White
- ISO size 9.0; from Orange to Yellow

Each department will need to consider locally how this change will be implemented, as it may take time for existing stock to be used.



*The four sizes of Intersurgical Guedel airways that have had their colours changed to conform to the updated standard.*

Product Code	Intersurgical Size	ISO Size	Previous Intersurgical Colour		New ISO Compliant Colour	
1000035 (S*)	000	3.5	Pink		Light Green	<b>NEW</b>
1100050 (S*)	00	5.0	Blue		Blue	
1110055 (S*)	0	5.5	Grey		Grey	
1111065 (S*)	1	6.5	White		Brown	<b>NEW</b>
1111570 (S*)	1.5	7.0	Yellow		White	<b>NEW</b>
1112080 (S*)	2	8.0	Green		Green	
1113090 (S*)	3	9.0	Orange		Yellow	<b>NEW</b>
1114100 (S*)	4	10.0	Red		Red	
1115120 (S*)	5	12.0	Purple		Purple	

This colour change effects all of Intersurgical's Guedel Airway products, including bulk-packed options.

(S\*) Add an S to the seven digit code number for the sterile version of this product eg. 5000000S

For more information, please visit [www.intersurgical.com/info/guedelairway](http://www.intersurgical.com/info/guedelairway)



Quality, innovation and choice

Interact with us



[www.intersurgical.com](http://www.intersurgical.com)