




# Our quick sizing guide for comparison

Patient	Age	Weight (Kg)	Traditional mask size	Quadralite masks	
Neonatal	Preterm	< 2	0		
	Term	2 – 4			
Infant	3 months	6	1		
	6 months	8			
	9 months	9			
Paediatric	1 year	10	2		
	2 years	12			
	3 years	14			
	4 years	16			
	5 years	18			
Small Adult	6 years	20	3		
	7 years	23			
	8 years	26			
	9 years	28			
	10 years	32			
	11 years	35			
	12 years	39			
Medium Adult	13 years	44	4		
	14 years	50			
	15 years	54			
	16 years	58			
Large Adult	16+ years	60 – 120	5		
Extra Large Adult	16 + years	> 120	6		