




## Nasal cannulae

A range of products for all patient groups



Oxygen and Aerosol Therapy • Nasal Cannulae

## Nasal cannula

### Conventional nasal cannula


A range of products for adult, paediatric and neonates are available with a wide selection of prong shapes to suit patient requirements. All oxygen tubing in the range is manufactured in a 'multi-channel' design, which permits oxygen flow even if the tubing is kinked.

### Also available in Satin


Cannula made from Intersurgical's soft satin material.

### Prong shape


The following prong shapes are available:




Curved prongs



Flared prongs



Curved and Flared prongs



Straight Prongs



Curved prongs with Intersurgical EarGuard™

### Intersurgical EarGuard™

The Intersurgical EarGuard option includes soft silicone foam cushions, which may reduce the risk of trauma to the patient's ears.

### Non-phthalate nasal cannula

All of the nasal cannula featured in this information sheet are now manufactured using non-phthalate materials to address DEHP exposure concerns.





### Multi-Channel oxygen tubing

Allows oxygen to continue to flow if the tubing becomes kinked

### Non-slip toggle

Enables user to adjust for best possible fit

### Non-phthalate

Addresses DEHP exposure concerns

### Clear material

Unobtrusive

### Over-ear design

Secure positioning

(Intersurgical EarGuard™ option available for added comfort)

### Soft prongs

Minimise irritation


### Location saddle

Assists in the optimum positioning of the prongs

Make an enquiry

Watch the video

Adult range of nasal cannulae



1165002

Code	Description	Length	Box Qty.
1161000	Adult, nasal cannula with straight prongs and tube	1.8 m	50
1161004	Adult, nasal cannula with straight prongs and tube	2.1 m	50
1161008	Adult, nasal cannula with straight prongs and tube	1.8 m	280
1162000	Adult, nasal cannula with straight prongs and tube	5 m	20
1162004	Adult, nasal cannula with straight prongs and tube	3.5 m	25
1165001	Adult, nasal cannula with curved prongs and tube	2.1 m	50
1169000	Adult, nasal cannula with straight prong, headset and tube		100
1165000	Adult, nasal cannula with curved prongs and tube	1.8 m	50
1161020	Adult, Satin nasal cannula with straight prongs and tube	1.8 m	50
1165020	Adult, Satin nasal cannula with curved prongs and tube	1.8 m	50
1165002	Intersurgical EarGuard™, adult, nasal cannula with curved prongs and tube	2.1 m	50
1165003	Intersurgical EarGuard, adult, nasal cannula with curved prongs and tube	5 m	40
1165021	Adult, Satin nasal cannula with curved prongs and tube	2.1 m	50
1167000	Adult, nasal cannula with curved/flared prongs and tube	1.8 m	50
1168000	Adult, nasal cannula with curved/flared prongs headset		100

Neonatal and paediatric range of nasal cannulae

Code	Description	Length	Box Qty.
1163000	Paediatric, nasal cannula with curved prongs and tube	2.1 m	50
1164000	Neonatal, nasal cannula with curved prongs and tube	2.1 m	50

IS10.8 • Issue 13 02.25



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK  
 T: +44 (0)118 965 6300 info@intersurgical.com www.intersurgical.com



The manufacturer Intersurgical Ltd is certified to ISO 14001:2015, ISO 9001:2015, ISO 13485:2016 and MDSAP

**Please think before you print**  
 Save energy and paper.  
 If you must print this information sheet please print it double sided.