

# Oral Care

the extended product range



## OroCare™ Aspire



OroCare™ Aspire is an innovative toothbrush with suction facility designed for use in intensive care units. It features a small brush head, equally suitable for adults and children. The soft bristles allow for proper oral hygiene and are designed to aid in the removal of plaque. The long neck allows the user to reach all areas of the oropharynx without interfering with ET tubes and other equipment.

OroCare™ Aspire has a built-in direct attachment to a suction line. This eliminates the need for an additional suction yankauer. Suction may conveniently be controlled by the centrally located thumb hole.



## OroCare™ Sensitive



OroCare™ Sensitive is a unique suction wand for cleaning delicate gums and tissues around the oropharynx. The integral moulded soft head provides a safer alternative to common suction swabs that can become detached. The soft atraumatic head may be used in circumstances where a toothbrush would be too irritating.

The design of the head allows application of antibacterial solutions, whilst it gently cleans the surface of the teeth and soft palate. A suction tube may be attached to the integrated connector and suction can be controlled by the conveniently located thumb hole.




## The OroCare™ suction line splitter



The OroCare™ suction line splitter has been designed for the simultaneous use of two suction lines on a single canister. This may become necessary when a closed suction system is already attached to the vacuum source.

# Ordering information

### Kits & toothbrushes

Product		Code	Box quantity
<b>OroCare™ Aspire suction toothbrush</b>		3016000	25
<b>OroCare™ Sensitive oral suction wand</b>		3017000	15
<b>OroCare™ suction line splitter</b>		3022000	25

IS1.37 • Issue 1 07.13